

Lancets

Clinical safety lancets

Lancing devices

Pen needles

Heel blood sampling devices

# stat<sup>®</sup> AUTO™ / Clinical Safety Lancets



**Tender**  
1.8mm depth / 26g.  
SAT-100

**Comfort**  
1.8mm depth / 21g.  
SAC-100

**Normal**  
2.4mm depth / 21g.  
SAN-100

**Xtra**  
2.8mm depth / 21g.  
SAX-100

**High-Flow**  
3.0mm depth / 21g.  
SAE-100

**Neonatal/Heelstick**  
1.2mm depth / 18g.  
SAH-100

- Simple one-step operation
- Non-reuse safety feature
- Home-use alternative for vision or dexterity-impaired
- Multi-language packaging
- Packaged 100 per box / 20 boxes per case
- Sterile
- Custom colors and packaging available
- US Patent #8,118,825



design  
manufacturing  
distribution

**U.S./International (001) 305 / 945 0011**  
**stat@statdevices.com**  
**statdevices.com**

STERILE | R



SMD-SAFS Rev. 6 6 / 16