

- Lancets
- Clinical safety lancets
- Lancing devices
- Pen needles**
- Heel blood sampling devices

stat® **SUPER-FINE™** / Insulin Pen Needles



enlarged
to show
detail

Mini	Lite	Micro	Comfort	Xtra	Normal	Standard
33g./4mm	32g./4mm	31g./5mm	31g./6mm	31g./8mm	30g./8mm	29g./12.7mm
PNI-100	PNL-100	PNM-100	PNC-100	PNX-100	PNN-100	PNS-100

- Silicone-coated, precision-engineered needle tips
- Each pen needle has inner sheath and outer cap
- Available in 7 sizes
- Multi-language packaging
- For use with most currently available insulin pens
- Packaged 100 per box - 96 boxes per case
- Sterile
- Private labeling available



design
manufacturing
distribution

U.S./International (001) 305 / 945 0011
stat@statdevices.com
statdevices.com

STERILE EO

